

A revolutionary way to reduce pain and swelling through the hot, cold, contrast and compression therapies you trust.

put a wrap on pain



VITALWEAR

VITALWRAP®
SYSTEM

Physicians and patients are fighting pain with precision, ease and flexibility.

Real relief, right where you need it. Reduce pain and swelling anywhere in the body-without drugs-in one breakthrough device. The **SADMERC** - approved VitalWrap System with its innovative family of wraps delivers targeted relief through hot, cold, contrast and compression therapies with precision and comfort.

SIMPLE:

Portable, easy-to-use design for patient self-treatment.

PRECISE:

Control temperature with one dial, 40° - 105° F.

FLEXIBLE:

Hot, cold, contrast and compression therapies.

COMPLETE:

Family of therapeutic wraps treats multiple body areas.



Wide-ranging applications — One versatile solution

Chronic Arthritis • Back pain • Sprains, strains and inflammation • Post-surgical recovery
Pain management • Functional mobility • Many other conditions



VITALWEAR



put a wrap on pain

Systems

order no.

00002 VitalWrap System 120 Vac, 60 Hz
(control unit, tubing set, standard VW)

00013 VitalWrap System 120 Vac, 60 Hz
(control unit, tubing set, without standard VW)

00123 VitalWrap System 220 Vac, 50 Hz
(control unit, tubing set, standard VW)

00165 VitalWrap System 220 Vac, 50 Hz
(control unit, tubing set, without standard VW)

Wraps

00255	Medium VitalWrap (4.25" x 52")	for elbow and calf
00001	Standard VitalWrap (6" x 60")	general purpose wrap for calf, knee, hamstring, hip, back, chest and abdomen
00118	Large VitalWrap (8" x 72")	broad coverage wrap for leg and back
00388	Back VitalWrap	for use on back
00119	Medium Shoulder VitalWrap	for use on either shoulder
00395	Large Shoulder VitalWrap	for use on either shoulder
00394	Wrist VitalWrap	for hand and wrist
00396	Hand VitalWrap	for arthritic hands
00397	Medium Knee VitalWrap	for post-surgery and CPM applications
00392	Full Knee VitalWrap	for post-surgery and CPM applications
00398	Ankle VitalWrap	for foot and ankle
00122	Cervical VitalWrap	for neck applications
00967	Hip II VitalWrap	for use on hip
01061	Cap VitalWrap	for use on head

VitalPad® made exclusively for XBACK® Solution ES/Compression

01434 Medium Vital Pad® for XBACK® Solution ES / Compression
Fits XBACK® Solution ES SM, MED | XBACK® Compression SM, MED

01435 Large Vital Pad® for XBACK® Solution ES / Compression
Fits XBACK® Solution ES LG, XL, 2XL | XBACK® Compression LG

01436 XLarge Vital Pad® for XBACK® Solution ES / Compression
Fits XBACK® Solution ES 3XL | XBACK® Compression XL, 2XL

01335 VitalPad for XBACK Solution, 4XL | for XBACK® Solution ES 4XL
XBACK® Bracing Services, Inc. is registered trademark.

Accessories

00094	Carrying Case	for easy portability of the VW System
00034	VitalWrap Tubing Set	
00987	VitalWrap "Y" Tubing Set	for use on 2 limbs
01417	IcePak™	reusable ice pack for cold therapy
00262	Hospital Grade Power Cord	

SADMERC Durable Medical Equipment Codes

code

E0217 Water circulating heat pad with pump

E0218 Water circulating cold pad with pump

E0249 Pad for water circulation heat unit

VitalWear, Inc.

384 Oyster Point Blvd., Suite 16
South San Francisco, CA 94080-1968

Tel: 1-800-553-4081
www.vitalwear.com

SYSTEM SPECIFICATIONS:

PHYSICAL	
size (approximate)	8" diameter x 13.5" height (20.3 cm x 34.3 cm)
weight (dry)	5 lbs (2.27 kg)
single tubing set	10.5 ft. (3.2m)
CONTROL SYSTEM	
type	manual user control
THERMAL SYSTEM	
range	40° - 105° F (4.4° - 40.6° C)
contrast therapy transition time	1 min typical from extremes in temperature range
operating period	2-6 hours typical continuous use
thermal cutout	113° F (45° C)
CIRCULATING SYSTEM	
reservoir capacity	1.1 gal (4.2 L)
reservoir fluid	ice water
flow rate thru VW	5.5 gph (21 Lph)
maximum pressure	18 psi
ELECTRICAL SYSTEM:	
voltage	120Vac, 60Hz (VIT-00002) or 220-240Vac, 50Hz, (VIT-00123)
internally fused	120 Vac, 60Hz, 300W F1, F2: 125V, 4A F3: 125V, 300µA 220-240Vac, 50Hz, 3.15A F1, F2: 250V, 2.5A F3: 250V, 300µA
leakage current	under 300 µA
OPERATING ENVIRONMENT:	
atmospheric pressure	525 to 795 mmHg
humidity	30% to 70% relative non-condensing
temperature	10°C to 40°C
TRANSPORT & STORAGE:	
atmospheric pressure	179 to 760 mmHg
humidity	10% to 95%
temperature	0°C and +50°C
REGULATORY:	
classification	class I equipment
type of equipment	type BF
REGULATORY APPROVALS:	
compliance	UL 60601-1, 1st edition, 2006-04-26: Standard for Medical Equipment Part 1: General Requirements for Safety. VITALWEAR™ and VITALWRAP™ are registered trademarks of VitalWear, Inc. Covered by US Patents: 7,211,104 7,191,798 7,658,205 7,694,693

For more information or to order VitalWrap Systems, wraps or accessories, contact our sales representatives at 1-800-553-4081 or visit www.vitalwear.com.